



# N-DEX® Plus # 8005PF

100% NITRILE GLOVE  
DISPOSABLE 8 MIL, 9.5" GLOVE  
NFPA 1999 COMPLIANT EMS GLOVE  
POWDERED OR POWDER-FREE



## Disposable Glove:

### Characteristics:

- 8005PF is Class One Medical Glove
- 8005 is EN Cat III Glove
- Sizes:

Small	8005(PF)S
Medium	8005(PF)M
Large	8005(PF)L
Extra-Large	8005(PF)XL

Cat. III EN 388



0123

2000

#### Performance Levels:

Abrasion: Level 0-4  
Cut Resistance: Level 0-5  
Tear: Level 0-4  
Puncture: Level 0-4

EN-374-3  
2003



J (Heptane)

0 1 2 3 4 5 6

K (Sodium Hydroxide 40%)

0 1 2 3 4 5 6

(Hydrochloric Acid 37%)

0 1 2 3 4 5 6

## ANSI/ISEA 105-2005 Performance

**Abrasion**

0 1 2 3 4 5 6

**Cut Resistance**

0 1 2 3 4 5

**Puncture**

0 1 2 3 4 5

## Description

Whether you're double-gloving over cut-resistant hand protection for meat processing applications or rebuilding the engine of a 1957 Chevy, striking the balance between proper protection and tactile sensitivity is critical. That's why you'll find the N-DEX® Plus on the hands of EMT, industrial and meat processing professionals with jobs that require the utmost in barrier protection, chemical and puncture resistance and tactile discrimination. The 8 mil disposable nitrile glove provides wearers barrier protection and comfort unequalled by supported nitrile gloves nearly twice its thickness. Available in lightly powdered and powder-free versions, N-DEX Plus special construction includes Best's low-modulus formulation, maximizing comfort and minimizing hand fatigue.

## Suggested Applications

- Emergency Medical Operations
- Manufacturing
- Food Processing
- Laboratories
- Janitorial
- Paint Applications
- Fuel Handling and Refineries

This list of suggested applications is partial and based on laboratory testing of glove performance. For a complete list of suggested applications, visit our website:

[www.showabestglove.com](http://www.showabestglove.com)

USA:  
Best Glove  
579 Edison Street  
Menlo, GA 30731  
Phone: 800-241-0323  
Fax: 888-393-2666  
Email: [USA@bestglove.com](mailto:USA@bestglove.com)

Canada:  
Best Glove Manufacturing, Ltd.  
253 Michaud Street  
Coaticook, Quebec J1A 1A9 Canada  
Phone: 800-565-BEST  
Fax: 819-849-6120  
Email: [Canada@bestglove.com](mailto:Canada@bestglove.com)

Europe:  
Best Manufacturing Europe N.V.  
Draaiboornstraat 6 bus 12  
2160 Wommelgem Belgium  
Phone: +32 (0)3 458 33 33  
Fax: +32 (0)3 458 34 57  
Email: [Europe@bestglove.com](mailto:Europe@bestglove.com)

Australasia:  
Best Manufacturing Australia Pty. Ltd.  
PO Box 187, Bexley  
NSW 2207 Australia  
Phone: +61 (0)2 8509 5401  
Fax: +61 (0)2 8509 5403  
Email: [Australiansales@bestglove.com](mailto:Australiansales@bestglove.com)